



MJ / GLOBAL SUPPLY

Precision for everyday practice.

Canadian surgical instruments catalog

REUSABLE • MANUAL • NON-POWERED

Focused essentials for clinical teams.

MJ Global Supply is developing a Canadian range of reusable manual surgical instruments for clinics, surgical centres, healthcare distributors and procurement teams. The collection begins with familiar general-use patterns and a controlled technical file for each approved product.

12

Class I* Canada

launch products

proposed category

Mississauga, Ontario

01

Controlled product selection

General-use manual patterns chosen for broad clinical familiarity.

02

Documented materials

Instrument-specific material evidence and lot traceability required before release.

03

Functional inspection

Alignment, ratchet, surface finish, dimension and cutting or gripping checks by lot.

04

Reprocessing information

Validated cleaning, disinfection and steam-sterilization instructions for reusable models.

Products 01-03

Common North American names and sizes. Final specifications remain subject to approved samples and manufacturer documentation.

SCISSORS MJG-SC-101 Operating Scissors Straight; sharp/sharp or sharp/blunt		SCISSORS MJG-SC-102 Mayo Dissecting Scissors Straight or curved; heavy pattern	
COMMON SIZE	5½ in / 14 cm	COMMON SIZE	6¾ in / 17 cm
PRIMARY USE	General cutting	PRIMARY USE	Firm tissue and sutures

SCISSORS MJG-SC-103 Metzenbaum Scissors Curved; blunt/blunt; delicate pattern	
COMMON SIZE	7 in / 18 cm
PRIMARY USE	Delicate dissection

Pre-launch product information

Products are not offered for sale until applicable licensing, final classification, bilingual labelling, specifications and technical documentation are complete.

Products 04-06

Common North American names and sizes. Final specifications remain subject to approved samples and manufacturer documentation.

HEMOSTASIS MJG-HF-201 Halsted Mosquito Forceps Straight or curved; serrated jaws		HEMOSTASIS MJG-HF-202 Kelly Hemostatic Forceps Straight or curved; transverse serrations	
COMMON SIZE	5 in / 13 cm	COMMON SIZE	5½ in / 14 cm
PRIMARY USE	Fine vessel clamping	PRIMARY USE	General hemostasis

FORCEPS MJG-TF-301 Dressing Forceps Straight; serrated; non-toothed	
COMMON SIZE	6 in / 15 cm
PRIMARY USE	Dressings and materials

Pre-launch product information

Products are not offered for sale until applicable licensing, final classification, bilingual labelling, specifications and technical documentation are complete.

Products 07–09

Common North American names and sizes. Final specifications remain subject to approved samples and manufacturer documentation.

FORCEPS MJG-TF-302 Adson Tissue Forceps Straight; 1 × 2 teeth COMMON SIZE 4¾ in / 12 cm PRIMARY USE Secure tissue handling	FORCEPS MJG-TF-303 Allis Tissue Forceps Ratchet handle; 4 × 5 teeth COMMON SIZE 6 in / 15 cm PRIMARY USE Firm tissue grasping
SUTURING MJG-NH-401 Mayo-Hegar Needle Holder Ratchet handle; serrated jaws COMMON SIZE 6 or 7 in PRIMARY USE General suturing	

Pre-launch product information

Products are not offered for sale until applicable licensing, final classification, bilingual labelling, specifications and technical documentation are complete.

Products 10–12

Common North American names and sizes. Final specifications remain subject to approved samples and manufacturer documentation.

SUTURING MJG-TC-501 Backhaus Towel Clamp Perforating tips; ratchet handle COMMON SIZE 5¼ in / 13 cm PRIMARY USE Drape fixation	RETRACTORS MJG-RT-601 Army-Navy Retractor Double-ended; handheld COMMON SIZE Standard pattern PRIMARY USE Shallow exposure
RETRACTORS MJG-RT-602 Weitlaner Retractor Self-retaining; 3 × 4 teeth COMMON SIZE 5½ in / 14 cm PRIMARY USE Maintained exposure	

Pre-launch product information

Products are not offered for sale until applicable licensing, final classification, bilingual labelling, specifications and technical documentation are complete.

Preparing the Canadian launch.

Buyer feedback, sample-evaluation discussions and distributor inquiries are welcome during product development.

SALES	sales@mjglobalsupply.ca
ORDERS	orders@mjglobalsupply.ca
REGULATORY	regulatory@mjglobalsupply.ca
ACCOUNTS	accounts@mjglobalsupply.ca

Activation status

These addresses are reserved in the launch plan and will activate after the business domain is purchased and connected. Replace this notice before public distribution.

*Classification note

Proposed classification only. Final classification follows the Canadian Medical Devices Regulations and each manufacturer's exact design and intended use. The cover photograph is a design presentation and does not represent exact production models.